

SHOULDER BRACE

QUICK START GUIDE

WHAT'S INCLUDED

- 1x Shoulder Brace
- 1x Elastic body strap
- 1x Short Strap Extender
- 1x Long Strap Extender

OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Vive Shoulder Brace as easy as possible. Check out the included link and QR code to help you through the process.



To see all of the FAQs in one place visit vhealth.link/d3cd9



HOW TO WEAR THE SHOULDER BRACE

1. Insert arm into brace.
2. Secure fastener strap around your chest.
3. Secure upper then lower strap around your bicep.

EXTEND THE LENGTH OF THE BODY STRAP

Your new shoulder support includes extenders for extending the length of the Body strap:

1. Secure both extenders to the corresponding hook and loop straps on the shoulder support.




2. Make sure the ends of the extenders are set so that they attach to the support. Adjust if necessary.
3. Put shoulder brace on and adjust as needed.

EXTEND THE LENGTH OF THE BICEP STRAP

Your new shoulder support includes 2 strap extenders for extending the size of the Bicep fastener:

1. Secure both extenders to the corresponding hook and loop straps on the shoulder support.



vive®

Distributed by

vive
health

8995 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com