

vive®



ICE THERAPY MACHINE

Owner's Manual
RHB1O49GRY

vivehealth.com

WHAT'S INCLUDED

- Ice Bucket with Lid
- Hose
- Universal Cold Therapy Pad
- 3x Adjustable Straps: Small, Medium and Large.
- Protective Barrier.



QUICK START GUIDE

Step 1:

Add ice to the bucket, filling to the ice fill line.



Step 2:

Add cold water to the bucket, filling to the water fill line.



Step 3:

Attach the lid to the top of the bucket.



Step 4:

Attach the hose to the lid, make sure there is an audible click.



Step 5:

Attach the cold therapy pad to the hose, make sure there is an audible click.



Step 6:

Plug in the power cord to the lid and an outlet. Then turn the unit on by flipping the switch next to the power plug.



Step 7:

Press and hold the "M" button to set the desired duration of cold therapy and the LPM (liter per minute). Use the up and down arrows to make changes and "M" to confirm. Press the "START" button.



COMMON APPLICATION AREAS

The cold therapy pad is suitable to be used on these areas.

1. Ankle



2. Calf



3. Knee



4. Quad



5. Hip



6. Shoulder



7. Elbow



8. Bicep



CARE AND MAINTENANCE

- Keep hoses and connectors free of dust and other contaminants.
- Store away from direct sunlight in a well ventilated area.

TROUBLESHOOTING

Reason	Solution
The pump will not turn on	Check all electrical connections and make sure the wall plug has power.
The cold therapy pad will not cool down	Make sure all trapped air is out of the pad before power is applied.
	Make sure the unit is completely filled with ice and water.
	Check all hose connections are secure and water is flowing.
	The pad is wrapped around your body part too tightly. Unwrap and rewrap with the unit running to ensure water is circulating.
The water is not flowing	Check the water level and add water if necessary.
	Check to make sure the hose is not kinked.

	<p>Make sure the connection between the hose and the pad is secure.</p>
	<p>The pad is wrapped around your body part too tightly. Unwrap and rewrap with the unit running to ensure water is circulating.</p>
<p>There is water leaking at the connector</p>	<p>Stop machine, disconnect the hose and reconnect listening for the "click". Restart the unit.</p>
	<p>Check the protective barrier used. If wet, replace with a new barrier.</p>
<p>There is water leaking at the cold therapy pad</p>	<p>The pad is defective. Replace immediately.</p>
	<p>Check the protective barrier used. If wet, replace with a new barrier.</p>
<p>There is magnetic or electromagnetic interference disrupting the unit</p>	<p>Move the unit away from the source of the interference.</p>
	<p>Change the power source away from the source of the interference.</p>

SPECIFICATIONS

- Dimensions: 11"x11"x11"
- Max water volume: 6L
- Unit Weight: 3.1kg (6.84lb)
- Timer settings: 5,10,15,20,25, 30 minutes
- Water flow -- 5 rates of flow (LPM)
- Water temperature indicator
- Celsius Range: 5 - 20 degrees
- Fahrenheit Range: 41 - 68 degrees
- 2 alarms (temperature and timer)
- Strap dimensions: 15", 27.75" and 46" in length
- Hose length: 5 foot length
- UL Listed
- FDA approved

QUESTIONS?

Our friendly customer service team is standing by to assist you. Feel free to connect on the phone, email or chat.



service@vivehealth.com



1-800-487-3808



vivehealth.com