Flexible & Versatile Alarm Settings

Patient Security Mode

Locks alarm limits, volume and time/date settings, prevents unintentional changes and provides automatic recall of locked limits.

Recall Function

Reinstatement of user's last customized alarm settings.

Veterinarian Use

7500FO could be used for veterinary applications requiring a pulse rate range of 300 bpm or less.

User-Defined Defaults

Allows user to adjust the alarm limit and volumedefault settings to their specific requirements.

Choose Your Options

Sensors

Use with NONIN's complete line of comfortable and reliable PureLight® sensors

Patient Cables

Optional use of right angle connector or sensor extension cables

Rolling Stand

The standard or deluxe 5-point stand increases portability

Carrying Case

Convenient for transport and storage

Fiberoptic Sensors

Use with the 7500FO only

Pole Mount System

Provides secure monitor placement

Data Management Software

7500 or 7500FO are compatible with nVISION®

Output Options

nVISION® Data Management Software for oximetry screening

PureSAT® Signal Processing Technology